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
The screenshot shows the LabSolutions main interface. At the top, there is a title bar with the text 'LabSolutionsメイン (System Administrator)' and a menu bar with 'ファイル(F)', '表示(V)', 'ウィンドウ(W)', and 'ヘルプ(H)'. Below the menu bar, there is a toolbar with a lock icon and a question mark icon. The main area is divided into a left sidebar and a main content area. The sidebar contains four buttons: '装置' (Devices), '解析ツール' (Analysis Tools), '管理ツール' (Management Tools), and 'マニュアル' (Manual). The '解析ツール' button is highlighted with a red border. The main content area is titled '解析ツール' and contains three icons: '再解析' (Re-analysis), 'ブラウザ' (Browser), and 'データベース' (Database). The '再解析' icon is highlighted with a red border.

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The screenshot displays a software interface for data analysis. On the left, a sidebar contains several menu items: 'メイン', 'データ解析', 'データレポート', 'メソッドに適用', 'ウィザード', and '手動波形処理' (highlighted with a red box). The main window is titled '再解析 (System Administrator) - [データ解析 - 学生実験 psi_b9_006.lcd]'. It features a file list on the left and a chromatogram plot on the right. The plot shows a signal with several peaks, with the most prominent one at approximately 160 minutes. The x-axis is labeled 'Time' and ranges from 0 to 375 minutes. The y-axis is labeled 'mV' and ranges from 0 to 12. The plot is titled 'クロマトグラムビュー' and includes a legend for '検出器A'. Below the plot, there is a 'メソッドビュー - 波形処理パラメータ' section with various settings for channel, width, slope, drift, and T. DBL. The bottom of the interface shows a 'メッセージ' window and a taskbar with the system clock at 15:09 on 2014/05/03.

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15:09
2014/05/03

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The screenshot shows a software application window titled '再解析 (System Administrator) - [データ解析 - 学生実験 psi b9_006.lcd]'. The 'ファイル(F)' menu is open, and the 'データレポート(I)' option is selected. A sub-menu is visible, showing '印刷(P) Ctrl+P' and 'PDF出力(O)...' as the active option. The main window displays two chromatograms. The top one is labeled 'クロマトグラムビュー' and shows a signal with several peaks. The bottom one is a zoomed-in view of the same signal, showing a large peak at approximately 15.5 minutes and two smaller peaks at approximately 8.5 and 19.5 minutes. The x-axis represents time in minutes (0.0 to 22.5) and the y-axis represents signal intensity in mV (0 to 12).

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The screenshot shows a software window titled "再解析 (System Administrator) - [データ解析 - 学生実験 psi b9_006.lcd]". The interface includes a menu bar (File, Edit, View, Method, Layout, Tools, Window, Help), a toolbar, and a sidebar with navigation options like "データレポート", "メソッドに適用", "ウィザード", and "手動波形処理". The main area displays a "クロマトグラムビュー" (Chromatogram View) with a plot of detector response (mV) over time (0 to 200 minutes). A "名前を付けて保存" (Save As) dialog box is open, showing the save location as "デスクトップ" (Desktop) and the file name as "学生実験 psi b9_006.pdf". The dialog also shows a list of system folders and a "保存する場所(D):" dropdown set to "デスクトップ". The file name field contains "学生実験 psi b9_006.pdf" and the file type is set to "PDF File (*.pdf)".

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